



MEL 90

The excimer laser that intelligently combines proven experience and exciting advancement



We make it visible.

The moment twice the speed puts
the finishing touch on proven experience.
This is the moment we work for.



// MEL 90
MADE BY ZEISS



MEL 90

My profile precisely

The use of excimer laser technology in refractive surgery starts in 1986 – with ZEISS. The MEL® 90 carries that legacy forward. Successfully combining years of expertise and innovative performance features, it offers users a highly rewarding handling experience.

Specially designed for the needs of the modern refractive surgeon, the MEL 90 incorporates proven safety as well as individually configurable options. Unique product qualities like the FLEXIQUENCE switch function, the new Triple-A profile and the outstanding intra-operative ablation speed of up to 1.3 seconds per diopter* create entirely new treatment prospects.

The MEL 90 from ZEISS – your preferences, your workflow, your expectations – precisely. The excimer laser that intelligently combines proven experience and exciting advancement.



** LASIK, myopia, 500 Hz, OZ: 6 mm*



MEL 90

Welcome

Dear Dr. Winkler

Arbeitsverordnungsstatus	
PKZ	200
USt	200
Gebühr	200

Welcome

Dear Dr. Winkler

Arbeitsverordnungsstatus	
PKZ	200
USt	200
Gebühr	200

My preferences precisely

When viewed closely, every refractive surgeon works somewhat differently, has personal preferences and needs regarding system handling and fosters his own surgical style. You will love the way the MEL® 90 takes all of these important factors into account.

Applied flexibility

FLEXIQUENCE

A true novelty – the FLEXIQUENCE switch function lets you flexibly choose between 250 Hz and 500 Hz frequency. It allows you to combine your previous clinical experience with excimer technology as well as the innovative advantages of the MEL 90. In fact, you can have twice the speed in a fraction of a second with a simple touch of the finger – enabling you to leverage both proven reliability and new treatment opportunities.

Predictable and gentle

Triple-A


The MEL 90 also features the Triple-A Advanced Ablation Algorithm. As a uniform ablation profile for a wide range of spherocylindrical (SCA) corrections, it simplifies treatment planning.

Triple-A offers a high degree of precision and predictability. In addition to intelligent energy correction, it comprises an aspherically optimized design that also focuses on minimal tissue removal. This results in gentle treatments of standard, but also eyes with higher and lower ametropias. All combined into one single ablation algorithm: Triple-A.



MEL 90 power

With its unique FLEXIQUENCE switch function (250 Hz and 500 Hz operation) and additional configuration options the MEL 90 is a true customized power package.



MEL 90 accuracy

The Triple-A profile enables fast, safe and highly precise treatments. Its excellent predictability, even with astigmatism and higher ametropias, offers additional safety.



MEL 90 time

When performing LASIK for myopia at 500 Hz, MEL 90 can intra-operatively ablate 1 diopter in as little as 1.3 seconds. Intuitive system guidance and speedy treatment planning allow for additional time savings.

My workflow precisely

You will especially value the treatment-oriented functions the MEL® 90 offers. Designed to make your work noticeably easier, they are also efficient and cost-effective. You can program the system to precisely fit your needs and that of your OR team – thereby adding focus to your overall workflow. You will love it.

Advanced system guidance

The systems' simple, intuitive graphic user interface supports fast treatment procedures. You can flexibly adjust the touch screen to wherever you need it and, thus, maintain a good ergonomic sitting posture throughout the operation.

Touch screen or keyboard control

For additional convenience, you can choose to make entries either per touch screen or keyboard.

Two workstations

Truly exceptional – a second, optional touch screen on the system, which allows the assistant to already prepare the next therapy while the surgeon is still concluding the post-operative care. Depending on the OR setup, the assistant touch screen can also be placed on a separate table.

Easy calibration

With the proven fluence test, laser calibration is very quick and easy.

Very service friendly

The MEL 90 is a highly reliable system requiring very little service maintenance. The unique beam path evacuation eliminates the need for flushing gas.

Seamless data integration

Patient data transfer is very easy thanks to the seamless interaction between the MEL 90 and other refractive platforms from ZEISS – VisuMax® and CRS-Master®.



MEL 90 connectivity

The workflow-optimized compatibility of the MEL 90, VisuMax and CRS-Master ensures a high level of practice consistency. An HD video port, network printer and PDF export capabilities offer additional workflow support.



MEL 90 teamwork

As a true team player, the MEL 90 greatly simplifies the interaction between the surgeon, assistant staff and system technology. When equipped with a second workstation, it can be configured to optimally complement the treatment routine of the OR team.



MEL 90 flexibility

Flexible touch screen positioning and optional keyboard control let you design your workstation just the way you want it.

My expectations precisely

Priority one, of course, is the satisfaction of your patients, which is largely dependent upon the safety and reliability of the treatment. The MEL® 90 combines high-quality elements and performance to help you optimally fulfill your customers' requirements. You and your patients will be delighted with the excellent clinical results.

Excellent precision

Triple-A corrections can be performed with impressive precision and a high degree reproducibility, even for very small and very large refractive errors, both with and without astigmatism. The new ablation profile and improved energy correction also offers other advantages such as few induced spherical aberrations and tissue-saving ablation.

High patient comfort

In combination, short treatments and seamless procedural flow help to enhance the patient experience. The comfortable patient supporting system ensures stress-free and stable positioning throughout the entire OR procedure.

Controlled OR environment

The ZEISS-exclusive sensor-controlled CCA+ plume removal system supports optimal atmospheric conditions throughout the treatment.

Safe eye tracking

The MEL 90 is equipped with an active eye tracker with an excellent response time to provide high-level treatment safety with very stable results.



MEL 90 intelligence

The Triple-A profile is based on an intelligent combination of refractive correction and clever energy distribution. It enables an aspherically optimized design with minimal ablation depth. As a result, treatments are faster and save tissue.

MEL 90 atmosphere

The unique, sensor-controlled CCA+ plume removal system supports optimal atmospheric conditions around the patient's eye – thereby also significantly improving the predictability of the results.

MEL 90 outcomes

The MEL 90 can be precisely tailored to the individual needs of the modern refractive surgical practice. Its preinstalled functions also offer surgeons optimal support for achieving excellent treatment results.

The refractive system landscape of ZEISS

As a comprehensive, well thought-out concept, ZEISS has put together an optimally matched refractive system landscape consisting of the MEL® 90 excimer laser, the VisuMax® femtosecond laser and the CRS-Master® treatment planning station. Operating in unison, these three leading edge products support a streamlined workflow and efficient patient management.

The MEL 90 performance spectrum

250 Hz

Femto-LASIK

LASIK

LASEK

PRK

PTK

500 Hz

Femto-LASIK

LASIK

- PRESBYOND Laser Blended Vision
- Individualized, topography-guided treatments
- Laser-independent sphero-cylindrical (SCA) treatment planning

VisuMax

Precision in all its facets

The VisuMax is a groundbreaking femtosecond system noted for its outstanding incision precision and speed as well as its gentle treatment qualities – whether for flap cuts, keratoplasty and incisions for ICR or the ZEISS-exclusive, minimally invasive method SMILE with ReLEx®.

CRS-Master

Simple treatment planning

The CRS-Master enables efficient treatment planning – even outside of the OR – for conventional as well as customized LASIK, Femto-LASIK, PRK and LASEK corrections. With its comprehensive capabilities, the CRS-Master is a valuable instrument for virtually any practice environment.



CRS-Master, VisuMax and MEL 90

The optimally interlinked refractive system landscape from planning to therapy

MEL 90, CRS-Master, PRESBYOND and ReLEx are not intended for sale in the United States.



MEL 90 flexibility for variable touch screen positioning

Integrated clipboard (optional) for patient documentation directly at the system

MEL 90 teamwork with the optional assistant workstation

Modern user interface for easy, intuitive guidance

Optional touch screen or keyboard operation

HD video port for high-resolution videos

Neat, tidy and safe cable routing

Additional storage compartment

MEL 90

ZEISS

PRESBYOND

Laser Blended Vision


The MEL 90 is also a first-rate platform for treating presbyopic patients. With the optional PRESBYOND® Laser Blended Vision module, ZEISS offers an advanced and true binocular laser correction method that far exceeds the limitations of conventional monovision.

The next level in eye care excellence

PRESBYOND Laser Blended Vision is an individualized treatment approach in which pre-operative spherical aberrations are used to calculate the ablation profile. The functional age of the eye can also be used for the calculation. This results in an increase in the depth of field of each individual eye, creating a customized fusion of the two images for near and distance vision – the so-called “Blend Zone.”

The method enables clear sight at all distances – near, far, even intermediate – and for a wide range of indications from -8.0 to +2.0 D. Emmetropic as well as astigmatic presbyopic patients (up to 2.0 cyl) can also be treated.

PRESBYOND Laser Blended Vision is tolerated by a higher number of patients than conventional monovision – up to 97%* as compared to only 59–67%** . In many cases, patients can read without glasses the very same day.



PRESBYOND Laser Blended Vision
Customized. All distances. Immediate.

* Reinstein DZ, Couch DG, Archer TJ. LASIK for Hyperopic Astigmatism and Presbyopia Using Micro-monovision With the Carl Zeiss Meditec MEL 80. *J Refract Surg.* 2009;25(1):37-58.

** Evans BJ. Monovision: a review. *Ophthalmic Physiol Opt.* 2007;27(5):417-439.



MEL 90

Technical data

Laser data

Type	ArF excimer laser (laser class 4 according to IEC 60825-1:2001 and IEC 60825-1:2007)
Wavelength	193 nm
Frequency	FLEXIQUENCE 250 Hz/500 Hz

Fixation laser

Type	Solid-state laser (laser class 1 according to IEC 60825-1:2001 and IEC 60825-1:2007)
Wavelength	532 nm

System data

Weight of MEL® 90	290 kg with gas cylinder
Weight of patient support system	232 kg
Dimensions of MEL 90 (W x D x H)	1630 mm x 730 mm x 1480 – 1700 mm
Power supply	100 V AC; 50/60 Hz; 17.5 A 120 V AC; 50/60 Hz; 14.6 A 208, 220, 230, 240 V AC; 50/60 Hz; 7.9 A
Approval	CE mark as per Medical Device Directive 93/42/EWG
Gas supply	Integrated ArF-Premix cylinder 10 l

Equipment

Surgical microscope	OPMI® pico with integrated HD video camera
Illumination	Ring illumination: stepless adjustment; sectional light; satellite illumination; optional slit lamp illumination
Active eye tracker	Infrared, pupil and limbus tracking, 1050 frames per second (fps), manual ablation center selection, automatic Pupil Center Shift Correction
CCA+ (plume removal system)	Integrated into the device, automatic adaptation for 250 Hz/500 Hz operation
Optional	Monitor with touch screen, keyboard, printer, CRS-Master, PRESBYOND® Laser Blended Vision

Spot scanning parameters

Beam dimensions	0.7 mm FWHM (full width at half maximum), Gaussian beam profile
-----------------	---

Phototherapeutic keratectomy

Area ablation	Programmable PTK shaping
---------------	--------------------------

Treatment range

According to CE guidelines	-12 D to +3 D (up to 3.0 D cyl)
Technically possible	-14 D to +8 D (up to 6.0 D cyl)



CE 0297



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.meditec.zeiss.com/contacts
www.meditec.zeiss.com/MEL90

Publication No: 000000-2054-914

EN_34_010_0041 Printed in Germany CZ-VIII/2013

The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. MEL, CRS-Master, PRESBYOND, VisuMax, ReLEx are either trademarks or registered trademarks of Carl Zeiss Meditec AG in Germany and/or other countries. Printed on elemental chlorine-free bleached paper.
© Carl Zeiss Meditec AG, 2013. All copyrights reserved.