



**VISALIS® 500**  
Accessory catalog



We make it visible.

# Table of contents

Page

I/A Tubing	3
<hr/>	
Cataract Accessories	3
<b>Phacoemulsification</b>	<b>3</b>
Handpiece	3
Phaco Sets	4
Tips	6
Additional Accessories	6
Sleeves	7
<hr/>	
<b>Irrigation/Aspiration</b>	<b>7</b>
Bimanual Technique	7
Coaxial Technique	8
<hr/>	
<b>Diathermy</b>	<b>10</b>
<b>Anterior Vitrectomy</b>	<b>10</b>
<hr/>	
Vitreoretinal Accessories	
<b>Trocar and Infusion Systems</b>	<b>11</b>
<b>Posterior Vitrectomy</b>	<b>11</b>
<b>Endoillumination</b>	<b>12</b>
<b>Injection and Removal Systems</b>	<b>13</b>
<b>Endodiathermy</b>	<b>14</b>
<hr/>	
Miscellaneous Accessories	14
System Options	15

## I/A Tubing

single-use		<p><b>I/A TUBING CASSETTE QUICKSET</b> 303060-0204-000 Cassette I/A tubing "QuickSet" with drainage bag (40 pcs)</p>
------------	---	--

## Cataract Accessories

### Phacoemulsification







#### Handpiece

reusable		<p><b>ULITE PHACO HANDPIECE</b> 303060-0001-000 Phaco Handpiece 4 piezoelectric, 40 KHz Tip stroke from 0 to 100 <math>\mu</math>m</p>
----------	---	--

# Cataract Accessories

## Phacoemulsification

### Phaco Sets







single-use		<b>PHACO SET 20G, U/S TIP 15° (10 pcs)</b>	20G Phaco Tip Straight, 15° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1000-000
single-use		<b>PHACO SET 20G, U/S TIP 30° (10 pcs)</b>	20G Phaco Tip Straight, 30° 20G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1001-000
single-use		<b>PHACO SET 20G, U/S TIP 30° FLARED (10 pcs)</b>	20G Phaco Tip Flared 30° 20G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.6-2.8	303060-1010-000
single-use		<b>PHACO SET 21G, U/S TIP 15° (10 pcs)</b>	21G Phaco Tip Biconical 15° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1002-000
single-use		<b>PHACO SET 21G, U/S TIP 30° (10 pcs)</b>	21G Phaco Tip Biconical 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1003-000
single-use		<b>PHACO SET 21G, U/S TIP 30° FLARED (10 pcs)</b>	21G Phaco Tip Flared 30° 21G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 2.2-2.5	303060-1011-000

\* Recommended incision size in mm

# Cataract Accessories

## Phacoemulsification

### Phaco Sets

single-use		<b>PHACO SET 22G, U/S TIP 15° (10 pcs)</b>  22G Phaco Tip Biconical 15° 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1004-000
single-use		<b>PHACO SET 22G, U/S TIP 30° (10 pcs)</b>  22G Phaco Tip 30° Biconical 22G Silicone Sleeve (2 pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1005-000
single-use		<b>PHACO SET 22G, U/S TIP 30° FLARED (10 pcs)</b>  22G Phaco Tip Flared 30° 22G Silicone Sleeve (2pcs) Test Chamber and Tip Wrench	RIS* 1.7-2.1	303060-1012-000
single-use		<b>20G SLEEVES, TEST CHAMBER SET (10 pcs)</b>  20G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.6-2.8	303060-1006-000
single-use		<b>21G SLEEVES, TEST CHAMBER SET (10 pcs)</b>  21G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 2.2-2.5	303060-1007-000
single-use		<b>22G SLEEVES, TEST CHAMBER SET (10 pcs)</b>  22G Silicone Sleeve (2 pcs) and Test Chamber	RIS* 1.7-2.1	303060-1008-000

\* Recommended incision size in mm



# Cataract Accessories

## Phacoemulsification

### Tips

reusable		19G PHACO TIP,15°	Straight 15°	303060-0050-000	RIS*	Sleeve 2.9-3.2	303060-0100-000
		19G PHACO TIP,30°	Straight 30°	303060-0051-000			
		19G PHACO TIP,45°	Straight 45°	303060-0052-000			
reusable		20G PHACO TIP,15°	Straight 15°	303060-0054-000	RIS*	Sleeve 2.6-2.8	303060-0101-000
		20G PHACO TIP,30°	Straight 30°	303060-0055-000			
		20G PHACO TIP,45°	Straight 45°	303060-0056-000			
reusable		20G PHACO TIP,BICONICAL	Biconical 15°	303060-0053-000	RIS*	Sleeve 2.6-2.8	303060-0101-000
reusable		20G PHACO TIP,20°,ANGLED	Angled 20°	303060-0062-000	RIS*	Sleeve 2.6-2.8	303060-0101-000
reusable		21G PHACO TIP,BICONICAL,15°	Biconical 15°	303060-0059-000	RIS*	Sleeve	303060-0102-000
		21G PHACO TIP,BICONICAL,30°	Biconical 30°	303060-0060-000			
reusable		22G PHACO TIP,BICONICAL,15°	Biconical 15°	303060-0063-000	RIS*	Sleeve	303060-0103-000
		22G PHACO TIP,BICONICAL,30°	Biconical 30°	303060-0065-000			
reusable		22G PHACO TIP,30°, SLEEVELESS BIMANUAL	Straight 30°	303060-0057-000	RIS*	1.0-1.2	
reusable		22G PROGR. PHACO TIP, SLEEVELESS BIMANUAL	Flared 30°	303060-0058-000	RIS*	1.0-1.2	
reusable		PARS PLANA PHACO TIP	20G (Posterior Fragmentation)	303060-0064-000			

### Additional Accessories

reusable		TEST CHAMBER (2 pcs)					303060-0202-000
reusable		PHACO TIP WRENCH					303060-0121-000

\* Recommended incision size in mm

# Cataract Accessories

## Phacoemulsification

### Sleeves

reusable		<b>19G SILICONE SLEEVE (GREY) 2.9-3.2MM</b> (3 pcs)	RIS* 2.9-3.2	303060-0100-000
reusable		<b>20G SILICONE SLEEVE (BLUE) 2.6-2.8MM</b> (3 pcs)	RIS* 2.6-2.8	303060-0101-000
reusable		<b>21G SILICONE SLEEVE (YELLOW) 2.2-2.5MM</b> (3 pcs)	RIS* 2.2-2.5	303060-0102-000
reusable		<b>22G SILICONE SLEEVE (GREEN) 1.7-2.1MM</b> (3 pcs)	RIS* 1.7-2.1	303060-0103-000

## Irrigation/Aspiration

### Bimanual Technique

single-use		<b>IRR &amp; ASP HANDPIECES BIMANUAL</b> (10 pcs) Irrigation Side Ports & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 2 Side Ports 0.45 mm		303060-0165-000
single-use		<b>IRR &amp; ASP HANDPIECES BIMANUAL FRONT</b> (10 pcs) Irrigation Front Port & Aspiration Bimanual Handpieces Aspiration: 21G, 0.35 mm Irrigation: 21G, 1 Front Port 0.5 mm		303060-0166-000
reusable		<b>IRRIGATION HANDPIECE, BIMANUAL</b> 21G, 0.6 mm		303060-0155-000
reusable		<b>ASPIRATION HANDPIECE, BIMANUAL</b> 21G, 0.3 mm		303060-0156-000

\* Recommended incision size in mm

## Cataract Accessories

### Irrigation/Aspiration

#### Coaxial Technique

single-use		<b>I/A HANDPIECE 20G STRAIGHT</b> 20G I/A Handpiece (10 pcs), Coaxial, Straight	RIS* 2.6-2.8	303060-0141-000
single-use		<b>I/A HANDPIECE 20G ANGLED 45°</b> 20G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.6-2.8	303060-0142-000
single-use		<b>I/A HANDPIECE 21G ANGLED 45°</b> 21G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 2.2-2.5	303060-0143-000
single-use		<b>I/A HANDPIECE 21G CURVED</b> 21G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 2.2-2.5	303060-0146-000
single-use		<b>I/A HANDPIECE 22G ANGLED 45°</b> 22G I/A Handpiece (10 pcs), Coaxial, Angled 45°	RIS* 1.7-2.1	303060-0144-000
single-use		<b>I/A HANDPIECE 22G CURVED</b> 22G I/A Handpiece (10 pcs), Coaxial, Curved 40°	RIS* 1.7-2.1	303060-0145-000
reusable		<b>I/A HANDPIECE, COAXIAL, REUSABLE</b> I/A Handpiece, Coaxial		303060-0140-000

\* Recommended incision size in mm



# Cataract Accessories

## Irrigation/Aspiration




### Coaxial Technique

reusable		1/A CANNULA, METAL SLEEVE, STRAIGHT, Ø1.6 MM		303060-0150-000
		1/A CANNULA, METAL SLEEVE, ANGLED, Ø1.6 MM		303060-0151-000
		1/A CANNULA, METAL SLEEVE, CURVED, Ø1.6 MM		303060-0152-000
		1/A CANNULA, MET. SLEEVE, 12 O'CLOCK, Ø1.6 MM		303060-0153-000
		1/A CANNULA, METAL SLEEVE, 90°, Ø1.6 MM		303060-0154-000
reusable		19G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.9-3.2	303060-0157-000
		19G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0158-000
reusable		20G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.6-2.8	303060-0159-000
		20G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0160-000
reusable		21G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 2.2-2.5	303060-0161-000
		21G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0162-000
reusable		22G I/A CANNULA, SILICONE SLEEVE, STRAIGHT	RIS* 1.7-2.1	303060-0163-000
		22G I/A CANNULA, SILICONE SLEEVE, ANGLED		303060-0164-000

\* Recommended incision size in mm

## Cataract Accessories

### Diathermy

reusable		<b>DIATHERMY PENCIL ERASER</b>	303060-0252-000
reusable		<b>DIATHERMY FORCEPS</b>	303060-0251-000
reusable		<b>DIATHERMY BIPOLAR CABLE</b>	303060-0250-000

### Anterior Vitrectomy




single-use		<b>20G A-VIT PROBE WITH SLEEVE (5 pcs)</b> 20G Anterior Vitrectomy Probe Silicone Sleeve	303060-0301-000
reusable		<b>20G A-VIT PROBE W/SLEEVE, REUSABLE</b> 20G Anterior Vitrectomy Probe Silicone Sleeve	303060-0300-000

## Vitreoretinal Accessories

### Trocar and Infusion Systems











single-use		<b>20G SCLERAL INFUSION CANNULA</b> (5 pcs)	303060-0503-000
single-use		<b>23G TROCAR KIT ONE-STEP</b> (10 pcs)	303060-0405-000
single-use		<b>23G TROCAR KIT ONE-STEP OIL INFUSION</b> (10 pcs)	303060-0406-000
single-use		<b>25G TROCAR KIT ONE-STEP</b> (10 pcs)	303060-0413-000
single-use		<b>25G TROCAR KIT ONE-STEP OIL INFUSION</b> (10 pcs)	303060-0412-000
single-use		<b>25G TROCAR KIT ONE-STEP</b> (5 pcs)	303060-0409-000

### Posterior Vitrectomy

single-use		<b>20G POSTERIOR VITRECTOMY PROBE</b> (5 pcs) High Speed	303060-0403-000
single-use		<b>23G POSTERIOR VITRECTOMY PROBE</b> (5 pcs) High Speed	303060-0407-000
single-use		<b>25G POSTERIOR VITRECTOMY PROBE</b> (5 pcs) High Speed	303060-0408-000

# Vitreoretinal Accessories

## Endoillumination

single-use		 <b>20G ENDO-ILLUMINATION PROBE (5 pcs)</b> 303060-0600-000
		 <b>20G WIDE-ANGLE ENDO-ILLUMINATION PROBE (5 pcs)</b> 303060-0601-000
single-use		 <b>20G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0602-000
		 <b>20G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs)</b> 303060-0603-000
single-use		 <b>23G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0604-000
		 <b>23G WIDE-ANGLE ENDO-ILLUMIN. MULTIPOINT (5 pcs)</b> 303060-0605-000
single-use		<b>25G ENDO-ILLUMINATION PROBE MULTIPOINT (5 pcs)</b> 303060-0606-000

# Vitreoretinal Accessories

## Injection and Removal Systems

single-use		<b>TWO WAY STOPCOCK (5 pcs)</b>	303060-0505-000
single-use		<b>ADAPTER AIR / CONTROLLED IRRIGATION (5 pcs)</b>	303060-0506-000
single-use		<b>AIR INJECTION TUBE WITH FILTER (5 pcs)</b>	303060-0508-000
single-use		<b>SILICONE-OIL INJECTION TUBE (5 pcs)</b>	303060-0509-000
single-use		<b>VISCOUS FLUID REMOVAL KIT (10 pcs)</b>	303060-0507-000

## Vitreoretinal Accessories


### Endodiathermy

single-use		20G ENDODIATHERMY PROBE (5 pcs)	303060-0262-000
single-use		23G ENDODIATHERMY PROBE (5 pcs)	303060-0263-000
single-use		25G ENDODIATHERMY PROBE (5 pcs)	303060-0264-000
reusable		DIATHERMY BIPOLAR CABLE	303060-0250-000





### Miscellaneous Accessories

single-use		SCREEN DRAPE VISALIS 500 (8 pcs)	303060-0352-000
single-use		REMOTE CONTROL COVER VISALIS (10 pcs)	303060-0353-000
single-use		TRAY COVER VISALIS 500 (2 pcs)	303060-0354-000

## Miscellaneous Accessories

reusable	 A rectangular blue and green sterilization tray with gold-colored contacts on the top surface.	<b>STERILIZATION TRAY</b> 26,5 x 16,3 x 4,6 cm	303060-0350-000
reusable	 A white, cylindrical touchscreen pen with a small circular button at the top.	<b>TOUCHSCREEN PEN (5 pcs)</b>	301640-0020-000

## System Options

	 A white, rectangular console tray with a metal handle on the right side and a small latch on the left.	<b>CONSOLE TRAY</b> with manually adjustable height	303041-9004-000
	 A white, rectangular video overlay box with a blue logo on the front panel.	<b>VIDEO OVERLAY BOX</b>	303041-9007-000
	 A white, rectangular IV-pole footswitch with two blue buttons and a cable attached to the right side.	<b>IV-POLE FOOTSWITCH</b>	303041-9006-000
	 A white, oval-shaped remote control with a cable attached to the top and a blue logo on the front.	<b>REMOTE CONTROL</b>	303041-9008-000

**Your local contact:****Argentina**

Carl Zeiss Argentina S.A.  
Calle Nahuel Huapi 4015 / 25  
C1430 BCO Buenos Aires  
Argentina  
Phone: +54 11 45 45 66 61  
bruzzi@zeiss.com.ar

**Australia**

Carl Zeiss Pty Ltd  
Tenancy Office 4, Level 1  
40-52 Talavera Road  
North Ryde NSW 2113  
Australia  
Phone: +61 2 9020 1333  
med@zeiss.com

**Austria**

Carl Zeiss GmbH  
Laxenburger Str. 2  
1100 Vienna  
Austria  
Phone: +43 1 79 51 80  
austria@zeiss.org

**Belgium**

Carl Zeiss NV-SA  
Ikaroslaan 49  
1930 Zaventem  
Belgium  
Phone: +32 2 719 39 11  
info@zeiss.be

**Brazil**

Carl Zeiss do Brasil Ltda.  
Av. Nações Unidas, 21711  
CEP04795-100 São Paulo  
Brazil  
Phone: +55 11 5693 5521  
medbrasil@zeiss.org

**Canada**

Carl Zeiss Canada Ltd.  
45 Valleybrook Drive  
Toronto, ON M3B 2S6  
Canada  
Phone: +1 800 387 8037  
czcmed@zeiss.com

**China**

Carl Zeiss Shanghai Co. Ltd.  
1/f., Ke Yuan Building  
11 Ri Yin Nan Road  
Waigaoqiao Free Trade Zone  
2005 Yang Gao Bei Road  
Shanghai 200131  
China  
Phone: +86 21 5048 17 17  
sro@zeiss.com.cn

**Czech Republic**

Carl Zeiss spol. s.r.o.  
Radlická 14/3201  
150 00 Prague 5  
Czech Republic  
Phone: +420 233 101 221  
zeiss@zeiss.cz

**France**

Carl Zeiss Meditec France SAS  
60, route de Sartrouville  
78230 Le Pecq  
France  
Phone: +33 1 34 80 21 00  
med@zeiss.fr

**Germany**

Carl Zeiss Meditec VG mbH  
Carl-Zeiss-Strasse 22  
73447 Oberkochen  
Germany  
Phone: +49 7364 20 6000  
vertrieb@meditec.zeiss.com  
Surgical Ophthalmology:  
Phone: +49 800 470 50 30  
iol.order@meditec.zeiss.com

**Hong Kong**

Carl Zeiss Far East Co. Ltd.  
Units 11-12, 25/F  
Tower 2, Ever Gain Plaza  
No. 88 Container Port Road  
Kwai Chung  
Hong Kong  
Phone: +852 2332 0402  
cffe@zeiss.com.hk

**India**

Carl Zeiss India Pvt. Ltd.  
Plot No.3, Jigani Link Road  
Bommasandra Industrial Area  
Bangalore - 560 099  
India  
Phone: +91 80 4343 8000  
info@zeiss.co.in

**Italy**

Carl Zeiss S.p.A.  
Viale delle Industrie 20  
20020 Arese (Milan)  
Italy  
Phone: +39 02 93773 1  
infomed@zeiss.it

**Japan**

Carl Zeiss Meditec Japan Co. Ltd.  
Shinjuku Ku  
Tokyo 160-0003  
22 Honchio-Cho  
Japan  
Ophthalmic instruments:  
Phone: +81 3 33 55 0331  
medsales@zeiss.co.jp  
Surgical instruments:  
Phone: +81 3 33 55 0341  
cmskoho@zeiss.co.jp

**Malaysia**

Carl Zeiss Sdn Bhd.  
Lot2, Jalan 243/51 A  
46100 Petaling Jaya  
Selangor Darul Ehsan  
Malaysia  
Phone: +60 3 7877 50 58  
malaysia@zeiss.com.sg

**Mexico**

Carl Zeiss de México S.A. de C.V.  
Avenida Miguel Angel de Quevedo  
496  
04010 Mexico City  
Mexico  
Phone: +52 55 59 99 0200  
cz-mexico@zeiss.org

**Netherlands**

Carl Zeiss B.V.  
Trapezium 300  
Postbus 310  
3364 DL Sliedrecht  
Netherlands  
Phone: +31 184 43 34 00  
info@zeiss.nl

**New Zealand**

Carl Zeiss (N.Z.) Ltd.  
15B Paramount Drive  
P.O. Box 121 - 1001  
Henderson, Auckland 0650  
New Zealand  
Phone: +64 9 838 5626  
med@zeiss.com

**Poland**

Carl Zeiss sp. z o.o.  
ul. Lopuzanska 32  
02-220 Warsaw  
Poland  
Phone: +48 22 858 2343  
medycyna@zeiss.pl

**Singapore**

Carl Zeiss Ptd. Ltd.  
50 Kaki Bukit Place  
Singapore 415926  
Singapore  
Phone: +65 6741 9600  
info@zeiss.com.sg

**South Africa**

Carl Zeiss (Pty.) Ltd.  
363 Oak Avenue  
Ferndale  
Randburg 2194  
South Africa  
Phone: +27 11 886 9510  
info@zeiss.co.za

**South Korea**

Carl Zeiss Co. Ltd.  
Seoul 121-828  
Mapo-gu  
141-1, Sangsu-dong  
2F, BR Elitel Bldg.  
South Korea  
Phone: +82 2 3140 2600  
korea@zeiss.co.kr

**Spain**

Carl Zeiss Meditec Iberia S.A.U.  
Ronda de Poniente, 15  
Tres Cantos  
28760 Madrid  
Spain  
Phone: +34 91 203 37 00  
info@zeiss.es

**Sweden**

Carl Zeiss AB  
Tegeluddsvaegen 76  
10254 Stockholm  
Sweden  
Phone: +46 84 59 25 00  
info@zeiss.se

**Switzerland**

Carl Zeiss AG  
Feldbachstrasse 81  
8714 Feldbach  
Switzerland  
Phone: +41 55 254 7534  
info@zeiss.ch

**Thailand**

Carl Zeiss Thailand  
90 CyberWorld Tower A,  
36<sup>th</sup> Floor, Unit A 3601  
230 Ratchadapisek Road  
Huaykhwang, Bangkok 10310  
Thailand  
Phone: +66 2 248 8787  
thailand@zeiss.com.sg

**United Kingdom**

Carl Zeiss Ltd.  
509 Coldhams lane  
Cambridge  
CAMBS  
CB1 3JS,  
United Kingdom  
Phone: +44 1707 871200  
info@zeiss.co.uk

**United States of America**

Carl Zeiss Meditec, Inc.  
5160 Hacienda Drive  
Dublin, CA 94568  
USA  
Phone: +1 925 557 4100  
info@meditec.zeiss.com

EN\_30\_010\_01131 Printed in Germany CZ-XI/2012  
The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information.  
Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.  
VISALUS is a registered trademark of Carl Zeiss. © 2011 by Carl Zeiss Meditec AG. All copyrights reserved.

